

Pressure Equalization Pads

with Tempur-Pedic[®] Technology





Optimized Patient Care and Surgical Support

Only Mizuho OSI® can deliver a comprehensive system of care to simplify your pressure management initiatives. From the surgical table to pressure redistributing surfaces and soft patient care kits, we deliver seamless solutions for optimal surgical positioning and patient support.

Mizuho OSI® **Pressure Equalization Pads** provide support, positioning and pressure redistribution during a surgical procedure. They are seamlessly integrated as a system of care with surgical tables and devices. We understand that quality materials, superior design and intimate knowledge of the surgical table translate into better patient outcomes.



Surgical Table
Manufacturer



Patient
Care Kits



Pressure
Equalization Pads




Pressure Equalization Pad Features

- Equalizes pressure load against the body
- Reduces incidence of hospital acquired pressure injuries (HAPI)¹
- Enhances patient comfort
- Adjusts firmness to body temperature
- Hypo-allergenic
- Latex-free
- Material is radiolucent, MRI safe
- Waterfall edges
- Moisture resistant
- Non-skid bottom
- Easy cleaning with standard disinfectants

Tempur-Pedic® vs. Conventional Foam

Studies show Tempur-Pedic® material provides up to 4 times the pressure reduction of conventional polyurethane in controlled testing.²

1. Moisture impervious, vapor permeable membrane
2. Flame resistant barrier mattress assembly
3. Tempur 110kg/cubic meter density for peak loading force distribution
4. Tempur-HD soft, dense absorption reduces localized pressure concentration



Medical Pads



Trios® and MTS

Part Number	Description	Size
Spinal Surgery Tops		
5939	Combination Hip Thigh Pad Set (2 each)	-
6946	Pillow Standard Size (Sling Pad)	-
6947	Pillow Queen Size (Sling Pad)	-
5379-7722	5379 Articulating/5376 Translating, Arm Board Pad	1"
5840-7030	Support Section (Leg Board Pad)	-
5840-7569	SST Single Chest Pad	-
Comfort Pads with Base		
6956	Comfort 8-Piece Pad Set	-
6956-3	Comfort Hip Pad, Large	-
6956-4	Comfort Hip Pad, Medium	-
6956-5	Comfort Hip Pad, Small	-
6956-6	Comfort Thigh Pad, Right	-
6956-7	Comfort Thigh Pad, Left	-
Comfort Pads without Base		
6956-31	Comfort Hip Pad without Base, Large	-
6956-41	Comfort Hip Pad without Base, Medium	-
6956-51	Comfort Hip Pad without Base, Small	-
6956-61	Comfort Thigh Pad without Base, Right	-
6956-71	Comfort Thigh Pad without Base, Left	-
Radiolucent Imaging Tops		
5840-7028	5927 and 7927 Table Pad	2"
Ortho Trauma Tops		
5855-7297	5855 Ortho Trauma Top Pad	-
5855-7357	Lower Leg Support Pad	-
5855-7367	Well Leg Pad	-
5855-7377	5855 Transfer Board Pad	-
5855-286	GM Board Right Non-Tempur Foam	-
5855-287	GM Board Left Non-Tempur Foam	-
5855-309	Lateral Pad Assembly Non-Temper Foam	-
Wilson Plus™		
5319-37	Wilson Plus Tempur-Pedic® Pads	-



Hana® / Hana SSXT® / ProFx®

Part Number	Description	Size
6875-7129	6875 Perineal Post Table Pad	-
6875-7120	6875 Hana Table Pad	-
6850-7600	6850 ProFx and Hana Arm Board Pad	-
6850-7590	6850-170 Leg Support Pad	-
6850-7560	6850 Perineal Post Table Pad	-
6850-7570	6850 ProFx Table Pad	-
6800-7140	6850-6000 ProFx Auxiliary Imaging Top/6800 and 6862 Allegro Pad, also fits Allegro 6800 and 6862	-
6850-7580	Lateral Perineal Support	-
6870-39	Hana SSXT Body Table Pad	-
6870-40	Hana SSXT Perineal Table Pad	-
5855-7377	Transfer Board Pad	-



ProAxis®

Part Number	Description	Size
6988-758	Chest Pad with Mounting Plate	-
5379-7722	Pad for 5379 and 6376 Arm Boards	-
6988-381	Hip/Thigh Pad with Mounting Bracket, Left	-
6988-389	Hip/Thigh Pad with Mounting Bracket, Right	-
5840-7030	Leg Support Pad	-
6988-442	Supine Top Pad	2"



Allegro®

Part Number	Description	Size
6800-7140	6800, 6862 Allegro / 6850-6000 ProFx Auxiliary Imaging Top, Table Pad	-
5356-7000	6900-12 Radiolucent Arm Board Pad	-
6900-25-7020	6900-25 AV Fistula Board Pad	-
6900-51-7001	6900-51 Universal Foot Board Pad	-

Positioning Pads

Positioning Pads with Tempur-Pedic® are available in six diverse shapes to meet your challenging positioning needs. These positioning devices quickly envelop contact points with a fluid density that reduces peak pressures, yet holds the patient in the desired position. With closed seams and no zipper these devices can be utilized in any position. Available individually or as a complete set.



Positioning Pads

Photo	Part Number	Description	Size
A - F	6950	6950-01 through 6950-06 OR Pad Set (6 pcs)	-
A	6950-01	Knee Block A Pad	10" x 8" x 8"
B	6950-02	Knee Block B Pad	14" x 14" x 4"
C	6950-03	Basic Roll Pad	18" x 6"
D	6950-04	Spinal Wedge Pad	24" x 10" x 6.75"
E	6950-05	Cradle Block Pad	20" x 16" x 6"
F	6950-06	Rectangular Pad	14" x 9" x 3"

Pillows & Arm Board Pads



Pillows

Part Number	Description	Size
6946	Pillow Standard Size (Sling Pad)	-
6947	Pillow Queen Size (Sling Pad)	-
6947-50	Disposable Pillow Case for 6946/6947 Pillow	-



Arm Board Pads

Part Number	Description	Size
5356-7000	5356 Arm Board Pad	24" x 6" x 2"
5366-7000	5356 Arm Board Pad	24" x 6" x 3"
5857-7000	5857 Cross Arm Support Pad	19" x 5" x 1"
5364-7018	Ultra Shoulder Translating Arm Board Pad	20" x 6" x 2"

GST Pads

All General Surgical Table Pads with Tempur-Pedic® are three piece sets with perineal cut-out. Individual pad sections can be ordered by placing an 'H' for Head, 'B' for Back or 'F' for Foot at the end of the pad set part number.

Due to manufacturer variables please confirm dimensions of existing pad(s) before ordering.



Steris® / AMSCO®

Part Number	Description	Size
Steris 3080 and 3085		
69402	3-Piece Pad with Perineal Cut-out	2"
69403	3-Piece Pad with Perineal Cut-out	3"
Steris 4085 and 5085		
CT-356079106	Steris 4085, 3-Piece Pad, 23.5" Foot Section with Cut-out	2"
CT-356079083	Steris 4085, 3-Piece Pad, 23.5" Foot Section with Cut-out	3"
CT-356079092	Steris 5085, 3-Piece Pad, 23" Foot Section with Cut-out	2"
CT-356079093	Steris 5085, 3-Piece Pad, 23" Foot Section with Cut-out	3"

MAQUET®

Part Number	Description	Size
CT-356129066	Maquet Alphamaxx, 3-Piece Pad Set with Perineal Cut-out	3"
CT-356129071	Maquet Alphastar Plus, 3-Piece Pad Set with Perineal Cut-out	2"

Skytron®

Part Number	Description	Size
6001, 6002, 6500 and 6600		
69412	3-Piece Pad with Perineal Cut-out	2"
69413	3-Piece Pad with Perineal Cut-out	3"
6700, 6701 and 6702		
69472	3-Piece Pad with Perineal Cut-out	2"
69473	3-Piece Pad with Perineal Cut-out	3"
3501, 3600B, 3601B and 3602		
69482	3-Piece Pad with Perineal Cut-out	2"
69483	3-Piece Pad with Perineal Cut-out	3"
3500 and Steris C-Max 4085 (No Foot Cutout)		
69492	3-Piece Pad with Perineal Cut-out	2"
69493	3-Piece Pad with Perineal Cut-out	3"

Stryker® / Berchtold®

Part Number	Description	Size
Berchtold B810		
69422	3-Piece Pad with Perineal Cut-out	2"
69423	3-Piece Pad with Perineal Cut-out	3"
694235	3-Piece Pad with Perineal Cut-out	3.5"
Operon D850, D860 and D820		
CT-356169011	D850, 4-Piece Pad Set with Perineal Cut-out	2"
CT-356169014	D850, 4-Piece Pad Set with Perineal Cut-out	3"

GST Arm Board Pads



Flat Arm Board Pads

Part Number	Description	Size
24"		
69441	Flat Arm Board Pad	24" x 6" x 1"
69442	Flat Arm Board Pad	24" x 6" x 2"
69443	Flat Arm Board Pad	24" x 6" x 3"
69444	Flat Arm Board Pad	24" x 6" x 4"
26"		
694461	Flat Arm Board Pad	26" x 6" x 1"
694462	Flat Arm Board Pad	26" x 6" x 2"
694463	Flat Arm Board Pad	26" x 6" x 3"
694464	Flat Arm Board Pad	26" x 6" x 4"



Contoured Arm Board Pads

Part Number	Description	Size
24"		
69452	Contoured Arm Board Pad	24" x 6" x 2"
69453	Contoured Arm Board Pad	24" x 6" x 3"
69454	Contoured Arm Board Pad	24" x 6" x 4"
26"		
694562	Contoured Arm Board Pad	26" x 6" x 2"
694563	Contoured Arm Board Pad	26" x 6" x 3"
694564	Contoured Arm Board Pad	26" x 6" x 4"

Miscellaneous



Wheel Chair Pads

Part Number	Description	Size
69461	Pressure Equalization Wheel Chair Pad	16" x 16" x 3"
69463	Pressure Equalization Wheel Chair Pad	18" x 18" x 3"
69464	Pressure Equalization Wheel Chair Pad	20" x 18" x 3"



Replacement Pads

Part Number	Description	Size
3369-1	Adult Perineal Post Pad, 6" Diameter	6"
5806-7010	5806 Radiolucent Imaging Top Pad	58" x 20" x 2"
5880-7150	Supine Operating Top Pad	40.5" x 20" x 2"
5847-218	Pelvic Reconstruction Kit, 6-Piece Pad Set	3"
CT-356149018	Stretcher Pad with Beveled Corners	76" x 31.5" x 3.5"

Montreal Lateral Positioner Pads

5875-7008	Head Section	31" x 20" x 2"
5875-7009	Center Section	12" x 18" x 1"
5875-7010	Foot Section	27" x 20" x 2"
5275-7030	2" Hand Surgery Table Pad	30" x 15" x 2"
5275-7031	3" Hand Surgery Table Pad	30" x 15" x 3"

Discontinued Table Replacement Pads

Part Number	Description	Size
6977-7000	Axis System Supine Imaging Top	71" x 20.5" x 2"



For more information or to order contact your account manager or call 1-800-777-4674.

References

1. DeFloor, T., DeBacquer, D., and Grypdonck, M. (2005). The Effect of Various Combinations of turning and Pressure Reducing Devices on the Incidence of Pressure Ulcers, International Journal of Nursing Studies, 42 (2005), 37-46.
2. Support Surfaces: Viscoelastic Expanded Urethane versus Conventional Foam. Eric Flam, PhD, PE and Loretta Raab, RPh, CCP, ECPM Jan/Feb 2005, pp 29 - 34.



30031 Ahern Avenue
Union City, CA 94587-1234 USA
Telephone: 510-429-1500
Toll Free: 800-777-4674
Fax: 510-429-8500
Outside USA: +1-510-429-1500
mizuhosi.com

Tempur-Pedic® is a registered trademark of Tempur-Pedic North America, LLC • All other trademarks mentioned are trademarks of their respective owners or holders.

Note: Mizuho OSI is constantly improving its products. All specifications are subject to change without notice. Mizuho OSI guarantees these products to be free of defects in material and workmanship with the exception of products which are misused, altered or damaged. All dimensions listed are approximate. Mizuho OSI is a Delaware Corporation.