



new

Hana[®] AATHA Sterile Drape Kit

Simplifies Draping For
Anterior Approach (AA)
Total Hip Arthroplasty
(THA)



Hana® AATHA Sterile Drape Kit

The **Hana Sterile Drape Kit** for Anterior Approach to Total Hip Arthroplasty as described by Dr. Joel Matta provides a sterile field with integrated component covers to match the unique geometry of the Hana Table.

- Simplifies draping for AATHA on the Hana
- Allows sterile Femoral Hook Support to be placed on the sterile Femoral Lift Extension
- Designed by Mizuho OSI, the manufacturer of the Hana Orthopedic Table
- Provides sterile protection over the Femur Lift
- Two U drapes allow for optimized surgical window for AATHA or revision surgery
- Can be used for patient's left or right operative site
- Clear lower U drape allows visualization of all Hana leg spar controls: gross and fine traction, internal and external leg rotation and spar lock handle
- Clear lower drape allows visualization of patient's legs and feet in traction boot



Drape Features

1. Complete upper body coverage and barrier from anesthesia
2. Integrated femur lift cover
3. Integrated femoral hook support drape
4. Upper and lower U drape frames surgical window*
5. Visualize all Hana controls and patient lower extremities*

*Only available with 6855SDK-6

Ordering Information

Part Number	Description	Bulk	Upper AUD Drape	Lower U Drape	Upper/Lower Drapes
6855AUD-6	Hana Anterior Hip Upper Sterile Drape (AUD)	6/cs	•		
6855SDK-6	Hana Sterile AATHA Drape Kit	6/cs	•	•	•

The Hana Anterior Hip Upper Sterile Drape (AUD) is compatible with U drapes with a 5" (12.7cm) or larger fenestration.
6855SDK-6 Kit Includes: Hana Anterior Hip Upper Sterile Drape, Lower U Drape and two clear Upper/Lower Drapes.



30031 Ahern Avenue
Union City, CA 94587-1234 USA
Telephone: 510-429-1500
Toll Free: 800-777-4674
Fax: 510-429-8500
Outside USA: +1-510-429-1500
mizuho.com
newhipnews.com

Founding Member



Patents Pending

Note: Mizuho OSI is constantly improving its products. All specifications are subject to change without notice. Mizuho OSI is a Delaware Corporation.

2022 © Mizuho OSI NW0996 REV A