



Patient Helper

Designed For
Total Convenience



Patient Helper

3

Points of easy adjustment

The Mizuho OSI® Patient Helper is designed to allow the patient to take additional control of their own mobility while the nurse does less physical moving of the patient. Features include convenient caregiver access to adjustments and the ability to lock the trapeze in place where it is most beneficial for the patient.

Convenient

- Easy one-person adjustments and quick set up
- Upright swivels to either side for ease of use
- Mounts to the head-end bracket of the bed
- Low profile fits easily through doorways and elevators for transport

Patient Benefits


- Assists patient with repositioning in bed, transferring to transport platforms and bed entry/egress
- Provides patient with independence and additional help with mobility

Care Giver Benefits

- Ergonomic design for assisting patients
- Upright swivels full 180 degrees to either side
- Upright quickly removes for patient access in emergency situations

Cost Effective

- Patient Helper Sleeves can be permanently attached to hospital beds
- Better utilization of assets: Upright and Trapeze can be moved from bed to bed



Ordering Information

REF

| | |
|--------|-------------------------------------|
| 3423 | Patient Helper with 3423-1 Sleeve |
| 3423-1 | Patient Helper Mounting Sleeve |
| 3423-2 | Patient Helper Upright and Trapeze |
| 3427 | Patient Helper with Adapter Bracket |
| 3427-1 | Patient Helper Mount Assembly |


Compatibility: 3423, 3423-1, and 3423-2 fits Secure II®, S3 and Epic II, VersaCare®, CareAssist®, Advanta/Advanta™ II, Advance® 2000, Century™, Basic Care, LTC, 840/870 Resident. 3427, 3427-1 and 3423-2 fits TotalCare®.



30031 Ahern Avenue
Union City, CA 94587-1234 USA
Telephone: 510-429-1500
Toll Free: 800-777-4674
Fax: 510-429-8500
Outside USA: +1-510-429-1500
mizuho.com
newhipnews.com

EC REP

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Note: Mizuho OSI is constantly improving its products. All specifications are subject to change without notice. Patient Helper is designed to fit several Hill-Rom and Stryker bed models. May require additional brackets provided by the bed manufacturer. Advance®, CareAssist®, TotalCare® and VersaCare® are registered trademarks of Hill-Rom. Advanta™ and Century™, are trademarks of Hill-Rom. Secure® is a registered trademark of Stryker Medical. • Mizuho OSI is a Delaware Corporation. Manufactured in the USA.